

Novel methods to detect microparticles *and improved analysis with flow cytometry*

Edwin van der Pol^{1,2}

Guus Sturk¹, Ton van Leeuwen^{2,3},
Rienk Nieuwland¹

July 23rd, 2011



¹Laboratory Experimental Clinical Chemistry; ²Biomedical engineering and Physics,
Academic Medical Center, Amsterdam, The Netherlands

³Biomedical Photonic Imaging, University of Twente, Enschede, The Netherlands



ISTH2011

Disclosures for Edwin van der Pol


In compliance with AMA requirements, ISTH makes the following disclosures to the session audience:

Research Support/P.I.	No relevant conflicts of interest to declare
Employee	No relevant conflicts of interest to declare
Consultant	No relevant conflicts of interest to declare
Major Stockholder	No relevant conflicts of interest to declare
Speakers Bureau	No relevant conflicts of interest to declare
Honoraria	No relevant conflicts of interest to declare
Scientific Advisory Board	No relevant conflicts of interest to declare

Presentation includes discussion of the following off-label use of a drug or medical device:

<N/A>

Introduction



- body fluids contain cell-derived vesicles
- clinically relevant information
- problem: vesicle detection

Objective

- explore the ability of novel and conventional methods to detect the *size* and *concentration* of vesicles in suspension


Methods

- standard population¹ of
 - vesicles prepared from human cell-free urine ($n=5$)
 - mixture of polystyrene beads
- analyzed by
 - Transmission Electron Microscopy
 - Nanoparticle Tracking Analysis
 - Resistive Pulse Sensing
 - Flow cytometry

Transmission Electron Microscopy

Philips CM-10

According to manufacturer



Vesicles

Beads


Nanoparticle Tracking Analysis



- determines *size* and *concentration* of vesicles in suspension
- recently successfully extended with *fluorescence detection*²

Nanoparticle Tracking Analysis


Nanosight NS500



Vesicles

Beads


Resistive Pulse Sensing



- determines *size* and *concentration* of vesicles in suspension

Resistive Pulse Sensing

iZon qNano




Vesicles



Beads

Flow cytometry and the refractive index

Polystyrene bead





$n = 1.61$




Silica bead



$n = 1.45$





Vesicle



$n_{\text{inside}} = 1.38 \pm 0.02$


$n_{\text{membrane}} = 1.48$



100 nm




Flow cytometer calibration




- absolute scattering power (mW) is calculated by Mie theory to relate detected scattering power (a.u.) to the diameter of beads


Flow cytometry calibration




Flow cytometry detection limit




Flow cytometry detection limit



Flow cytometry detects *multiple vesicles as single count*




vesicles from human urine
filtered with 220 nm filter



89 nm silica beads at
concentration $1E10 \text{ ml}^{-1}$

Beam volume \approx
54 μ l



At a concentration
of $1E10 \text{ ml}^{-1}$, >800
vesicles are
simultaneously
present in the beam

Conclusion

- Nanoparticle Tracking Analysis and Resistive Pulse Sensing are promising methods to determine size and concentration of *single* vesicles in suspension (P-MO-405)
- detection of vesicles by flow cytometry is attributed to scattering from *multiple* vesicles (P-MO-404)